



make your patients' dreams come true



USP
ULTRA SHORT PULSE

..... todos los fototipos, todo el año



primelase®


Feel the new laser generation

 **COCOON**
MEDICAL



primelase[®]

The unique laser on the market with combined wavelengths



ULTRA SHORT PULSE

primelase®

New technology for hair removal

- Extreme fast and versatility
- Highest Safety
- No pain
- Best Results
- Maximum Profitability



primelase®

New technology for hair removal

All wavelenghts in
one system




USP
ULTRA SHORT PULSE

755nm

For treating skin phototypes I, II, III and IV
with blond and fine hair.

Phototypes



primelase[®]

New technology for hair removal

All wavelengths in
one system




USP
ULTRA SHORT PULSE



810nm

Recommended for treating all skin
phototypes.

Phototypes



primelase®

New technology for hair removal

All wavelengths in
one system




USP
ULTRA SHORT PULSE

Combined 810nm, 940nm and 1060nm



Recommended for III, IV, V and VI.

Phototypes



primelase®

Perfect combination

Versatility, effectiveness and safety




USP
ULTRA SHORT PULSE

1060nm

Indicated for treating dark phototypes IV, V and VI.

Phototypes




primelase[®]

New technology for hair removal

All wavelengths in one system

**Unique laser on the
market** with combined
wavelengths and high peak
power 4000W



primelase[®]

Versatility,
effectiveness
and safety

Spots size

				
				
Spot 10x9	Spot 10x9	Spot 30x17	Spot 20x9	Spot 11x9
Potencia de pico: 1800W	Potencia de pico: 2400W	Potencia de pico: 4800W	Potencia de pico: 4000W	Potencia de pico: 4000W
Spot 20x9	Spot 20x9			
Potencia de pico: 3600W	Potencia de pico: 4800W			

primelase®

Versatility

2 Works Mode Without Pain

Dynamic Mode
(5Hz and 10Hz)

pain free mode



primelase[®]

Versatility

2 Works Mode
Without Pain

Static Mode
(1Hz, 2Hz and 3Hz)




primelase®

Safety and Results

Advantages of perfect square pulse


- Best heating hair and skin not
- More effective
- Without Laser peaks (control laser radiation in the skin)
- Less exposure Laser
- More safety without side effects

primelase 10J/5ms



1° ms - 2J
2° ms - 2J
3° ms - 2J
4° ms - 2J
5° ms - 2J

others 10J/5ms




1° ms - 0,8J
2° ms - 1J
3° ms - 2,2J
4° ms - 3J
5° ms - 3J

primelase[®]

Highest Safety

Uniform in all treatments

- Uniform heating of the spot
- Uniform delivery of the energy in all cm² of the spot
- Warranty of homogeneous results
- Without side effect



primelase®

Safety-Whitout Pain Crystal Freeze

High skin protection

Always guarantees the same temperature +5°C


- Regardless of the room temperature
- Regardless of treatment speed
- Regardless of treatment fluence



The skin temperature is reduced by around 25 °C with the parallel cooling system *crystal freeze*

The skin increases to approximately 45 °C

The hair follicle reaches 70-80 °C generating a micro-burn at its base.



primelase[®]

Safety, Speed (Profitability)

Advantages for a big spots

Big Spot Size
=
Less Joules
needed

Short Pulse
Duration
=
Less Joules
needed

Combination of Big Spot Size + Short Pulse
=
Less Joules needed
=
Highest Safety with Maximum Results

- 2 month with out growing hair
- Lose 20%-30% hair per session

primelase[®]

Best Results
primelase effects



Treatment men
breast

Hair explosion




primelase[®]

Best Results primelase effects



Treatment men legs

Immediately after
treatment



primelase[®]

Best Results primelase effects



Treatment men
abdomen

Immediately after
treatment




primelase®

Best Results primelase effects




Treatment woman
abdomen

Immediately after
treatment



primelase[®]

Best Results
primelase effects




Treatment men stomach

Skin type VI

8 ms /16 Jcm²/ 8Hz /1060 nm

primelase[®]



primelase®

Speed and Profitability

The fastest on the market

Area	Time (min)
Legs	20
Arms	10
Back	5
Armpit	3



primelase®

Cost-Profitability

Why primelase?

- No consumables
- No maintenance
- Without fibre
- Without laser lamp
- Long life hand pieces
- Without special electrical supply needs



primelase[®]

Conclusion



Extreme
Versatility



Highest
Safety



Better
Effectiveness



Maximum
Profitability



Best Results



www.cocoonmedical.com

Follow us on:

